Tutoriel sur dcm4che2

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Version initiale 2009

Version finale 2016

Table des matières

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- ***org.dcm4che2.audit.log4j.helpers*** :

Class Summary:

SyslogWriter : SyslogWriter is a wrapper around the java.net.DatagramSocket class so that it behaves like a java.io.Writer.

- ***org.dcm4che2.audit.message*** :

Provides classes for generating XML formatted Audit Messages compliant to the "Audit Trail and Node Authentication (ATNA) Integration Profile" specified in the IHE IT Infrastructure Technical Framework, which itself utilizes RFC 3881 - Security Audit and Access Accountability Message XML Data Definitions for Healthcare Applications and its extension by DICOM Supplement 95: Audit Trail Messages.

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- **org.dcm4che2.imageioimpl.plugins.dcm**

Class Summary	
DicomImageReader	
DicomImageReaderSpi	
DicomImageWriter	Write DICOM files containing images - handles transcoding to a new representation, as well as transcoding the images.
DicomImageWriterSpi	Provide information about the dicom image writers

- **org.dcm4che2.io**

Class Summary	
ContentHandlerAdapter	
DataOutputStreamAdapter	This class adapts a regular output stream onto a data output stream, such as would be found in an ImageOutputStream
DicomInputStream	
DicomOutputStream	
ImageInputStreamAdapter	
SAXReader	
SAXWriter	
StopTagInputHandler	
TranscoderInputHandler	

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Class Summary	
Composite	
CRImage	
DXImage	The Digital X-Ray (DX) Image Information Object Definition specifies an image that has been created by a digital projection radiography imaging device.
Image	
SpatialFiducials	The Fiducials IOD specifies the spatial relationship between the Composite Fiducial instance, to one or more images.

- **org.dcm4che2.iod.module.composite**

This package contains classes which represent the Modules that are common to all Composite Image IODs (C.7 COMMON COMPOSITE IMAGE IOD MODULES).

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Class LookupTable		
Retour	Méthode	Explication
LookupTable	combine(LookupTable other, int outBits, Boolean inverse, short[] pval2out)	Création de la ColorModel
LookupTable	combine(LookupTable vlut, LookupTable plut, int outBits, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, non-linear Modality LUT and non-linear VOILUT.
LookupTable	createLut(int inBits, Boolean signed, int outBits, DicomObject mLut, DicomObject voiLut, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, non-linear Modality LUT and Window Center/Width.
LookupTable	createLut(int inBits, Boolean signed, int outBits, DicomObject mLut, DicomObject voiLut, DicomObject pLut, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, Rescale Slope/Intercept, non-linear VOILUT and non-linear Presentation LUT.
LookupTable	createLut(int inBits, Boolean signed, int outBits, DicomObject mLut, float center, float width, java.lang.String vlutFct, boolean inverse, short[] pval2out)	Create LUT for given i/o range, non-linear Modality LUT and Window Center/Width
LookupTable	createLut(int inBits, Boolean signed, int outBits, DicomObject mLut, float center, float width, java.lang.String vlutFct, DicomObject pLut, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, non-linear Modality LUT, Window Center/Width and non-linear Presentation LUT
LookupTable	createLut(int in Bits, Boolean signed, int outBits, float slope, float intercept, DicomObject voiLut, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, Rescale Slope/Intercept and non-linear VOI LUT.
LookupTable	CreateLut(int inBits, Boolean signed, int outBits, float slope, float intercept, DicomObject voiLut, DicomObject pLut, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, Rescale Slope/Intercept, non-linear VOI LUT and non-linear Presentation LUT

Class LookupTable		
Retour	Méthode	Explication
LookupTable	createLut (int inBits, Boolean signed, int outBits, float slope, float intercept, float center, float width, java.lang.String vlutFct, Boolean inverse, short[] pval2out)	Create ramp or sigmoid LUT for given i/o range, Rescale Slope/Intercept and Window Center/Width.
LookupTable	createLut (int inBits, Boolean signed, int outBits, float slope, float intercept, float center, float width, java.lang.String vlutFct, Boolean inverse, short[] pval2out)	Create LUT for given i/o range, Rescale Slope/Intercept, Window Center/Width and non-linear Presentation LUT.
LookupTable	createLutForImage (DicomObject img, DicomObject mlutObj, DicomObject voiObj, DicomObject pLut, float center, float width, java.lang.String vlutFct, short[] pval2out)	Create LUT for given DICOM image and output range.
LookupTable	createLutForImage (DicomObject img, int outBits, short[] pval2out)	Create LUT for given DICOM image and output range.
LookupTable	createLutForImageWithPR (DicomObject img, DicomObject mLut, DicomObject voiLut, boolean inverse, int outBits, short[] pval2out)	Create LUT for given DICOM image, non-linear VOI LUT and output range.
LookupTable	createLutFromWL (DicomObject img, DicomObject mLut, float center, float width, java.lang.String vlutFct, boolean inverse, int outBits, short[] pval2out)	Create LUT for given DICOM image, modality LUT, and Window Center/Width and output range.
int	getOffset	@return off
int	length()	@ return la taille d'un tableau
lookupTable	scale(int outBits, Boolean inverse, short[] pval2out)	
Int	toIndex(int in)	.

Class OverlayUtils		
Retour	Méthode	Explication
int	<u>extractFrameNumber</u> (int imageIndex)	Extra the frame number portion of the overlay number/imageIndex value.
Java.awt.image.BufferedImage	<u>extractOverlay</u> (DicomObject ds, int overlayNumber, javax.imageio.ImageReader reader, java.lang.String rgbs)	Read an overlay image or region instead of a regular image.
int	<u>getOverlayHeight</u> (DicomObject ds, int overlayNumber)	Reads the height of the overlay - needs to be done separately from the primary width even though they are supposed to be identical, as a stand alone overlay won't have any width/height except in the overlay tags.
boolean	<u>getOverlayWidth</u> (DicomObject ds, int overlayNumber)	Reads the width of the overlay - needs to be done separately from the primary width even though they are supposed to be identical, as a stand alone overlay won't have any width/height except in the overlay tags.
boolean	<u>isOverlay</u> (int imageIndex)	Returns true if the given frame number references an overlay - that is, is the form 0x60xx yyyy where xx is the overlay number, and yyyy is the overlay frame number.

Class VOIUtils		
Retour	Méthode	Explication
Int[]	calcMinMax(DicomObject img, java.awt.image.Raster raster)	Gets the min/max value from a data buffer, in the raw pixel data.
boolean	containsVOIAttributes(DicomObject dobj)	Return true if the specified object contains some type of VOI attributes at the current level
Float[]	getMinMaxWindowCenterWidth(DicomObject img, DicomObject pr, int frame, java.awt.image.Raster raster)	This method returns the minimum and maximum window center widths, based on the Modality LUT and image information.
boolean	isModalityLUTcontainsPixelIntensityRelationshipLUT(DicomObject img)	Return true if the specified image object contains a pixel intensity relationship LUT, based on SOP class
boolean	isModalityLUTcontainsPixelIntensityRelationshipLUT(java.lang.String uid)	Return true if the specified image object contains a pixel intensity relationship LUT
DicomObject	selectModalityLUTObject(DicomObject img, DicomObject pr, int frame)	Finds the applicable DicomObject containing the Modality LUT information for given frame within the image.
DicomObject	selectVoItemFromPr(java.lang.String uid, DicomObject pr, int Frame)	Searches for a Softcopy VOILUT Sequence for the given frame, or one for just the SOP instance
DicomObject	selectVoObject(DicomObject img, DicomObject pr, int frame)	Finds the application DicomObject containing the VOI LUT information for the given frame within the image

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Package org.dcm4che2.imageio.plugins.dcm

<i>DicomImageReadParam</i>	
DicomStreamMetaData	

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Package org.dcm4che2.imageioimpl.plugins.dcm

<i>DicomImageReader</i>	An abstract superclass for parsing and decoding of images.
DicomImageReaderSpi	Constructs a blank DicomImageReaderSpi.
DicomImageWriter	Write DICOM files containing images - handles transcoding to a new representation, as well as transcoding the images.
DicomImageWriterSpi	Provide information about the dicom image writers

Class DicomImageReader		
boolean	canReadRaster()	Returns true if this plug-in supports reading just a Raster of pixel data.
void	copyReadParam(javax.imageio.ImageReadParam src, javax.imageio.ImageReadParam dst)	
javax.imageio.ImageTypeSpecifier	createImageSpecifier()	Create an image type specifier for the entire image
void	dispose()	Allows any resources held by this object to be released.
javax.imageio.ImageReadParam	getDefaultReadParam()	
int	getHeight(int imageIndex)	Returns the height in pixels of the given image within the input source.
javax.imageio.metadata.IIOMetadata	getImageMetadata(int imageIndex)	Gets any image specific meta data.
java.util.Iterator<javax.imageio.ImageTypeSpecifier>	getImageTypes(int imageIndex)	Returns an Iterator containing possible image types to which the given image may be decoded, in the form of ImageTypeSpecifiers (dicom).
int	getNumImage(boolean allowSearch)	Returns the number of regular images in the study.
javax.imageio.metadata.IIOMetadata	getStreamMetadata()	Return a DicomStreamMetadata object that includes the DICOM header.

Class DicomImageReader(suite)		
Retour	Méthode	Explication
int	getWidth (int imageIndex)	Returns the width in pixels of the given image within the input source.
void	initImageReader (int imageIndex)	Sets the input for the image reader.
void	postDecompress ()	
Java.awt.image.BufferedImage	read (int imageIndex, javax.imageio.ImageReadParam param)	Reads the provided image as a buffered image.
byte[]	readBytes (int imageIndex, javax.imageio.ImageReadParam param)	Reads the bytes for the given image as raw data.
Java.awt.image.Raster	readRaster (int imageIndex, javax.imageio.ImageReadParam param)	Read the raw raster data from the image, without any LUTs being applied.
void	reset ()	Restores the ImageReader to its initial state.
void	setInput (java.lang.Object input, boolean seekForwardOnly, boolean ignoreMetadata)	Sets the input source to use to the given ImageInputStream or other Object.
int	getNumImages (boolean allowSearch)	Returns the number of regular images in the study.

Class DicomImageReaderSpi		
Retour	Méthode	Explication
boolean	canDecodeInput (java.lang.Object input)	Returns true if the supplied source object appears to be of the format supported by this reader.

javax.imageio.ImageReader	createReaderInstance (java.lang.Object extension)	Returns an instance of the ImageReader implementation associated with this service provider.
java.lang.String	getDescription (java.util.Locale locale)	Returns a brief, human-readable description of this service provider and its associated implementation.

Class DicomImageReaderSpi(suite)		
Retour	Méthode	Explication
boolean	canWriteSequence()	Writing as a sequence is actually the preferred approach.
javax.imageio.metadata.IIOMetadata	convertImageMetadata (javax.imageio.metadata.IIOMetadata metadata, javax.imageio.ImageTypeSpecifier type, javax.imageio.ImageWriteParam param)	The image metadata provided MUST be of the correct type for the child image writer, as the meta-data will only be written to the child image.
javax.imageio.metadata.IIOMetadata	convertStreamMetadata (javax.imageio.metadata.IIOMetadata metadata, javax.imageio.ImageWriteParam param)	Can't convert anything except existing DicomStreamMetadata
void	endWriteSequence()	Finish writing the header data to the stream.
byte[]	extractImageEncoding (javax.imageio.IIOImage iioimage, javax.imageio.ImageWriteParam param)	This method gets the encoded image from the given object as a byte array of data.
javax.imageio.metadata.IIOMetadata	getDefaultImageMetadata (javax.imageio.ImageTypeSpecifier image, javax.imageio.ImageWriteParam writeParam)	No easy way to figure out what this one is as the stream meta-data isn't yet available.
javax.imageio.metadata.IIOMetadata	getDefaultStreamMetadata (javax.imageio.ImageWriteParam arg0)	Get a default set of DICOM data to use in the stream meta-data.
void	prepareWriteSequence (javax.imageio.metadata.IIOMetadata metadata)	Start writing the initial sequence.

void	setupWriter (javax.imageio.metadata.IOMetadata metadata)	Sets up the child writer if it hasn't already been setup
void	updateDicomHeader (DicomStreamMetaData metadata, java.awt.image.BufferedImage image)	This updates the metadata with relevant information from the image, such as the size etc.
void	write (avax.imageio.metadata.IOMetadata metadata, javax.imageio.IIOImage iioimage, javax.imageio.ImageWriteParam param)	Writes the image, including the DICOM metadata.
void	writeBytesToSequence (_byte[] data, javax.imageio.ImageWriteParam param)	Write the given image as a byte array to the sequence.
void	writeToSequence (javax.imageio.IIOImage iioimage, javax.imageio.ImageWriteParam param)	Write the given image to the sequence.

Class DicomImageWriterSpi		
Retour	Méthode	Explication
boolean	canEncodeImage (javax.imageio.ImageTypeSpecifier type)	Indicate the DICOM can encapsulate the given type
javax.imageio.ImageWriter	createWriterInstance (java.lang.Object extension)	Create a dicom image writer
Java.lang.String	getDescription (java.util.Locale locale)	Get the description of this image writer type.

Package org.dcm4che2.imageioimpl.plugins.rle

RLEImageReader	
RLEImageReaderSpi	

Class RLEImageReader		
Retour	Méthode	Explication
int	getHeight (int imageIndex)	Returns the height in pixels of the given image within the input source.
javax.imageio.metadata.IIOMetadata	getImageMetadata (int imageIndex)	Returns an IIOMetadata object containing metadata associated with the given image, or null if the reader does not support reading metadata, is set to ignore metadata, or if no metadata is available.
java.util.Iterator<javax.imageio.ImageTypeSpecifier>	getImageTypes (int imageIndex)	Returns an Iterator containing possible image types to which the given image may be decoded, in the form of ImageTypeSpecifiers.
int	getNumImages (boolean allowSearch)	Returns the number of images, not including thumbnails, available from the current input source.
javax.imageio.metadata.IIOMetadata	getStreamMetadata (int imageIndex)	Returns an IIOMetadata object representing the metadata associated with the input source as a whole (i.e., not associated with any particular image), or null if the reader does not support reading metadata, is set to ignore metadata, or if no metadata is available.
int	getWidth (int imageIndex)	Returns the width in pixels of the given image within the input source.
java.awt.image.BufferedImage	read (int imageIndex, javax.imageio.ImageReadParam param)	Reads the image indexed by imageIndex and returns it as a complete BufferedImage, using a supplied ImageReadParam.
void	setInput (java.lang.Object input, boolean seekForwardOnly, boolean ignoreMetadata)	Sets the input source to use to the given ImageInputStream or other Object.

Class RLEImageReaderSpi

Retour	Méthode	Explication
boolean	canEncodeImage (javax.imageio.ImageTypeSpecifier type)	Indicate the DICOM can encapsulate the given type
javax.imageio.ImageWriter	createWriterInstance (java.lang.Object extension)	Create a dicom image writer
Java.lang.String	getDescription (java.util.Locale locale)	Get the description of this image writer type.

Package org.dcm4che2.io

ContentHandlerAdapter	Receive notification of the logical content of a document.
DataOutputStreamAdapter	This class adapts a regular output stream onto a data output stream, such as would be found in an ImageOutputStream
DicomInputStream	
DicomOutputStream	
ImageInputStreamAdapter	
SAXReader	
SAXWriter	
StopTagInputHandler	
TranscoderInputHandler	

Class ContentHandlerAdapter		
Retour	Méthode	Explication
void	characters (char[] ch, int start, int length)	Receive notification of character data dicom.
void	endElement (java.lang.String namespaceURI, java.lang.String localName, java.lang.String qName)	Receive notification of the end of an element
void	setDocumentLocator (org.xml.sax.Locator locator)	Receive an object for locating the origin of SAX document events.
Class DataOutputStreamAdapter		
Retour	Méthode	Explication
	org.xml.sax.Attributes atts)	

void	writes (byte[] b, int off, int len)	Writes len bytes from array b, in order, to the output stream.
void	write (int b)	Writes to the output stream the eight low-order bits of the argument b.

Class DicomInputStream		
Retour	Méthode	Explication
DicomObject	getDicomObject ()	Obtient un Object du dicom donc lecture des métadata
long	getEndOfFileMetaInfoPosition ()	Obtient la fin des métadatas du dicom
byte[]	getPreamble ()	
long	getStreamPosition ()	
TransferSyntax	getTransferSyntax ()	
int	level ()	
void	mark (int readlimit)	Marks the current position in this input stream.
int	read (byte[] b, int off, int len)	Reads up to len bytes of data from this input stream into an array of bytes. This method blocks until some input is available.
DicomObject	readDicomObject ()	Lecture des métadonnées dicom
void	readDicomObject (DicomObject dest, int len)	
DicomObject	readFileMetaInformation ()	Read File Meta Information from this stream. If there is no File Meta Information on current stream position, null is returned. Otherwise the stream will be parsed until the

		end of the File Meta Information is detected and a DicomObject containing the File Meta Information is returned.
void	readFully (byte[] b)	
void	readFully (byte[] b, int off, int len)	
int	readHeader ()	Lecture de l'entête.
void	readItem (DicomObject dest)	Lecture des Items
void	readItems (DicomElement sq, int sqLen)	
boolean	readValue (DicomInputStream dis)	Called by the input stream when reading a DICOM value in the stream.
void	reset ()	Repositions this stream to the position at the time the mark method was last called on this input stream.
void	setEndOfFileMetaInfoPosition (long fmi EndPos)	
void	setHandler (DicomInputHandler handler)	
void	setStreamPosition (long pos)	
long	skip (long n)	
DicomElement	sq ()	
int	tag ()	
long	tagPosition ()	

int	<code>valueLength()</code>	
VR	<code>vr()</code>	

Class DicomOutputStream		
Retour	Méthode	Explication
void	<code>close()</code>	Closes this file output stream and releases any system resources associated with this stream.
void	<code>finish()</code>	Finishes writing compressed data to the output stream without closing the underlying stream. Use this method when applying multiple filters, and the transfer syntax is a deflator transfer syntax.
byte[]	<code>getPreamble()</code>	.
long	<code>getStreamPosition()</code>	
TransferSyntax	<code>getTransferSyntax()</code>	
boolean	<code>isAutoFinish()</code>	Indicate if the stream is finished automatically (compressed data written) when the dataset is written
boolean	<code>isExplicitItemLength()</code>	
boolean	<code>isExplicitItemLengthIfZero</code>	
boolean	<code>isExplicitSequenceLength()</code>	

Class DicomOutputStream(suite)

Retour	Méthode	Explication
boolean	isExplicitSequenceLengthIfZero()	.
void	serializeDicomObject (DicomObject attrs)	Only for internal use by DicomObjectSerializer.
void	setAutoFinish (boolean autoFinish)	.Set to false to not auto finish the stream on writing a data set - useful for writing a DataObject followed by some additional DICOM that is custom written, eg images or related large data.
void	setExplicitItemLength (boolean explicitItemLength)	
void	setExplicitItemLengthIfZero (boolean explicitItemLengthIfZero)	
void	setIncludeGroupLength (boolean includeGroupLength)	Indicate if the stream is finished automatically (compressed data written) when the dataset is written
void	setPreamble (byte[] preamble)	
void	setStreamPosition (long pos)	
void	setTransferSyntax (java.lang.String tsuid)	
void	write (byte[] b, int off, int len)	Writes len bytes from the specified byte array starting at offset off to this file output stream.
void	write (int b)	Writes the specified byte to this output stream
void	writeCommand (DicomObject attrs)	
void	writeDataset (DicomObject attrs, java.lang.String tsuid)	Write a DICOM dataset to the output stream.

void	writeDataset (DicomObject attrs, TransferSyntax transferSyntax)	Write a DICOM dataset to the output stream.
void	writeDicomFile (DicomObject attrs)	Write a DICOM object to the output stream using the specified DicomObject to obtain the transfer syntax UID and other attributes.
void	writeFileMetaInformation (DicomObject attrs)	
void	writeHeader (int tag, VR vr, int len)	
void	writeItem (DicomObject item, TransferSyntax transferSyntax)	Write an item (DicomObject) to the output stream.

Class ImageInputStream		
Retour	Méthode	Explication
void	mark (int readlimit)	Marks the current position in this input stream.
boolean	markSupported ()	Tests if this input stream supports the mark and reset methods.
int	read ()	Reads the next byte of data from the input stream.
int	read (byte[] b,int off, int len)	Reads up to len bytes of data from the input stream into an array of bytes.
void	reset ()	Repositions this stream to the position at the time the mark method was last called on this input stream.
long	skip (long n)	Skips over and discards n bytes of data from this input stream.

Class SAXReader		
Retour	Méthode	Explication
DicomObject	readDicomObject()	Lit les métadatas du dicom
DicomObject	readDicomObject(DicomObject dcmObj)	

Class SAXWriter		
Retour	Méthode	Explication
java.io.File	getBaseDir()	
Int[]	getExclude()	
boolean	readValue (DicomInputStream in)	Called by the input stream when reading a DICOM value in the stream.
void	setBaseDir (java.io.File baseDir)	
void	setExclude (int[] exclude)	
void	write (DicomObject attrs)	

Class StopTagInputHandler		
Retour	Méthode	Explication
boolean	readValue (DicomInputStream in)	Called by the input stream when reading a DICOM value in the stream.

Class TranscoderInputHandler

Retour	Méthode	Explication
boolean	<code>readValue(DicomInputStream in)</code>	Called by the input stream when reading a DICOM value in the stream.

Package org.dcm4che2.data

This package contains lower level classes that deal with the DICOM Data Dictionary (PS3.6) Data Structure and Encoding (PS3.5).

Please note that some of these classes (such as [Tag](#)) are automatically generated by XSLT from DICOM Standard Text MS Word documents. So your amendments and/or modifications will be lost with the next run of dcm4che2/dcm4che2-core/dcm4che2-core-dict/build.xml to include attributes defined in final text supplements and correction items additional to Part 6 and 7 of the base standard text.

BasicDicomObject	
DateRange	
DicomObjects	
DicomObjectToStringParam	
ElementDictionary	
Implementation	
PersonName	
RessourceLocator	
SpecificCharacterSet	
Tag	Provides tag constants
TransferSyntax	
UID	Provides tag constants
UIDDictionary	
VR	This class provides enumeration of DICOM Value Representation
VRMap	

Class BasicDicomObject		
Retour	Méthode	Explication
boolean	accept (DicomObject.Visitor visitor)	<p>Calls @link{Visitor#visit} for each element in this Dataset. Returns false, if @link{Visitor#visit} returns false for any element.</p> <p>Parameters:</p> <p>visitor - <i>Visitor</i> object, which method @link{Visitor#visit} is called for each element in this Dataset.</p>
void	add (DicomElement a)	
boolean	bigEndian ()	Called by the input stream when reading a DICOM value in the stream.
void	bigEndian (boolean bigEndian)	
boolean	cacheGet ()	
void	cacheGet (boolean cacheGet)	
void	clear ()	Removes all elements from this Dicom Object.
DicomObject	command ()	
Java.util.Iterator<DicomElement>	commandIterator ()	Returns an iterator over Command elements (0000, eeee) in this Dicom Object.
boolean	contains (int tag)	Returns true, if this Dicom Object contains the specified Element.
boolean	containsAll (DicomObject keys)	Returns true if this Dicom Object contains all of the elements in the specified Dicom Object.
boolean	containsValue (int tag)	Returns true, if this Dicom Object contains the specified Element with a value length > 0.
void	copyTo (DicomObject dest)	

DicomObject	dataset()	
Java.util.Iterator<DicomElement>	datasetIterator()	Returns an iterator over Data elements (group tag > 2) in this Dicom Object.
boolean	equals(java.lang.Object o)	

Class BasicDicomObject(suite)		
Retour	Méthode	Explication
DicomObject	exclude(int[] tags)	
DicomObject	excludePrivate()	
DicomObject	fileMetaInfo()	
Java.util.Iterator<DicomElement>	fileMetaInfoIterator()	Returns an iterator over File Meta Information elements (0002, eeee) in this Dicom Object.
DicomElement	get(int tag)	
DicomElement	get(int[] tagPath)	
DicomElement	get(int[] tagPath, VR vr)	
byte[]	getBytes(int tag)	
byte[]	getBytes(int[] tagPath)	
byte[]	getBytes(int[] tagPath, boolean bigEndian)	
byte[]	getBytes(int tag, boolean bigEndian)	
java.util.Date	getDate(int tag)	
java.util.Date	getDate(int[] tagPath)	
java.util.Date	getDate(int[] tagPath, java.util.Date defVal)	

java.util.Date	getDate (int[] itemPath, int daTag, int tmTag)	
java.util.Date	getDate (int[] itemPath, int daTag, int tmTag, java.util.Date defVal)	
java.util.Date	getDate (int tag, VR vr)	
java.util.Date	getDate (int tag, VR vr, java.util.Date defVal)	
DateRange	getDateRange (int tag)	
DateRange	getDateRange (int[] tagPath)	
DateRange	getDateRange (int[] tagPath, DateRange defVal)	
DateRange	getDateRange (int[] tagPath, int daTag, int tmTag)	
DateRange	getDateRange (int[] tagPath, int daTag, int tmTag, DateRange defVal)	
DateRange	getDateRange (int[] tagPath, VR vr)	
DateRange	getDateRange (int[] tagPath, VR vr, DateRange defVal)	

Class BasicDicomObject(suite)		
Retour	Méthode	Explication
DateRange	getDateRange (int tag, DateRange defVal)	
DateRange	getDateRange (int daTag, int tmTag)	
DateRange	getDateRange (int daTag, int tmTag, DateRange defVal)	
java.util.date[]	getDates (int tag)	
java.util.date[]	getDates (int[] tagPath, java.util.Date[] defVal)	
java.util.date[]	getDates (int[] tagPath, VR vr, java.util.Date[] defVal)	
java.util.date[]	get (int[] tagPath, VR vr)	

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The number of bits used to encode the VR and VL fields depends on the transfer syntax and the VR. In the default implicit VR little endian transfer syntax, the VR is never sent and the VL is always 32 bits. In the all other transfer syntaxes, the VR is always sent and is 16 bits. For OB, OW, SQ, UT and UN VRs, the 16 bit VR is followed by 16 bits of zeroes, then a 32 bit VL. For all other VRs, a 16 bit VR field is followed by a 16 bit VL. The idea is to stay aligned on a 32 bit boundary when 32 bit VLs are required for "big" values, and save 32 bits for the others. In retrospect, this "optimization" has caused nothing but trouble, particularly when adding new VRs to the standard, or trying to send "big" values in "small" VRs.

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```
File sourceFile = new File("cardio.dcm");
```

```

ImageInputStream in =
ImageIO.createImageInputStream(sourceFile); //On crée une entrée de bit

```

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```
Iterator<ImageReader> iter = imageIO.getImageReadersByFormatName("DICOM");  
// returns an Iterator containing all currently registered ImageReaders that claim to be able to decode  
the named format.
```

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Signification:

Reads the provided image as a buffered image. It is possible to read image overlays by providing the 0x60000000 number associated with the overlay. Otherwise, the imageIndex must be in the range 0...numberOfFrames-1, or 0 for a single frame image. Overlays can be read from PR objects or other types of objects in addition to image objects. param can be used to sepecific GSPS to apply to the image, or to override the default window level values, or to return the raw image.

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```
/**  
 * Reads the provided image as a buffered image. It is possible to read  
 * image overlays by providing the 0x60000000 number associated with the  
 * overlay. Otherwise, the imageIndex must be in the range  
 * 0..numberOfFrames-1, or 0 for a single frame image. Overlays can be read  
 * from PR objects or other types of objects in addition to image objects.  
 * param can be used to sepecific GSPS to apply to the image, or to override  
 * the default window level values, or to return the raw image.  
 */
```

```
@Override  
public BufferedImage read(int imageIndex, ImageReadParam param)  
    throws IOException {  
    if (OverlayUtils.isOverlay(imageIndex)) {  
        readMetaData();  
        String rgbs = (param != null) ? ((DicomImageReadParam) param)  
            .getOverlayRGB() : null;  
        return OverlayUtils.extractOverlay(ds, imageIndex, this, rgbs);  
    }  
  
    initImageReader(imageIndex);  
    if (param == null) {  
        param = getDefaultReadParam();  
    }  
  
    BufferedImage bi;
```

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```
float[] emboss = new float[] { -2,0,0, 0,1,0, 0,0,2 };
```

- The Kernel class defines a matrix that describes how a specified pixel and its surrounding pixels affect the value computed for the pixel's position in the output image of a filtering operation. The X origin and Y origin indicate the kernel matrix element that corresponds to the pixel position for which an output value is being computed.

```
Kernel kernel = new Kernel(3, 3, emboss);
```

- Constructs a ConvolveOp given a Kernel.

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